



## UNI-132D Slide Warming Table



### Smart Waterbath for Today's Smart Labs

The Step Down slide warmer's unique design saves space and allows easy retrieval of slides since the end of the slide is not on the flat surface. Great for paraffin tissue section mounting in the fields of cytology, histology, pathology and biology. It is also an ideal choice for laboratories, education institutes and scientific research organisations. The heated surface leaves the label part of the slide protruding and therefore cool. The slide supports are angled to achieve a tiered lay-out, and the back legs of the instrument are higher to bring slides to the horizontal plane, which keeps the cover slips in place. The result is a space-saving slide warmer with any slide accessible and cool to pick-up.



## Features:

- The Step Down slide warmer is a "tiered" laboratory slide warmer that is easy to assemble and use.
- The Step Down has 5 levels holding up to 50 slides.
- The new design of the Step Down surface allows for easier and safer pick-up of slides.
- The slide warmer has an electrophoretic surface coating.
- The Step Down warms slides via heating elements which transfer the heat across the entire surface.
- The Step Down slide warmer is designed for use in the fields of cytology, histology, and biology.
- Independent Temperature Control.
- Temperature selection between Ambient to 90°C

## Product Specifications:

Dimensions: W x L x H (cm)  
403 x 520 x 150

Voltage:  
220V, 50 Hz, Single Phase, A.C.



UNIMEDITREK PRIVATE LIMITED  
14 Chitra Vihar,  
New Delhi-110092  
India  
Ph.: +91-11-65651688  
Mob.: +91-9810533445  
E-mail: sales@unimeditrek.com  
Website: www.unilabequipments.com